

Chromatin Immunoprecipitation Assay (ChIP Assay)**Day 1****Prep:**

- ChIP Dilution Buffer (CDB) (@ 4°C)
- liquid nitrogen in container
- rotator @ 4°C
- tubes (1 for blocked beads, 1 for each ChIP sample)
- heat block/water bath @37°C
- BSA (NEB) (@-20°C)
- Protease Inhibitor (PI) (@-20°C)
- Protein G agarose beads (@4°C)
- Antibody

Solutions:

1.) ChIP Dilution Buffer (CDB) (16.7 mM Tris-HCl pH 8.1, 167 mM NaCl, 0.01% SDS, 1.1% Triton-X 100)

1M Tris-HCl, pH 8.1	835 ul	
5M NaCl	1.67 ml	
10% SDS	50 ul	Keep on ice at all times.
10% Triton-X 100	5.5 ml	
H ₂ O	41.9 ml	
Total	50 ml	

2.) Protease Inhibitor (PI) (50X)

- Dissolve 1 tablet in 1000 ul ddH₂O.

Directions:**Preclear Chromatin**

- 1.) Gently shake and slightly warm agarose beads @ rm. temp. until well resuspended.
- 2.) In tube on ice, add: 1 ml CDB + 80 ul agarose beads slurry.
- 3.) Quickly thaw chromatin (@-80°C) in 37°C heatblock/waterbath. Place on ice.
- 4.) Add at least 5 ug chromatin to CDB and agarose beads slurry.
- 5.) Rotate @ 4°C for 30-60 minutes.

Block Agarose Beads (do at same time as Preclearing Chromatin)

- 1.) In tube, add: 1 ml CDB + 40 ul agarose beads slurry.
 - For each additional ChIP, add another 40 ul agarose beads slurry.
- 2.) Rotate @4°C for 5 minutes.
- 3.) Spin @ 3.4 x 1000 rpm (1.1 rcf) for 30 sec.
- 4.) Aspirate supernatant.
 - Be careful to not remove beads!
- 5.) Add 1 ml CDB.
- 6.) Repeat steps 2-5 for 2 more washes, but after 3rd wash, do not add CDB.
 - During 2nd wash step, thaw PI and BSA @ rm. temp.

7.) After last aspiration, add:

	<u>1X</u>	
(agarose beads only)	(20 ul)	(washed)
CDB	68 ul	
10 mg/ml BSA	10 ul	
<u>50X PI</u>	<u>2 ul</u>	
Total	100 ul	

8.) Rotate blocked beads overnight @ 4°C.

Antibody Binding

1.) After 30-60 min. rotation @ 4°C, on ice, carefully remove precleared chromatin supernatant (~1 ml) from beads and put into new 1.5 microcentrifuge tube.

2.) Add:

- 20 ul PI (50X)
- ~2 ug antibody

3.) Rotate overnight @ 4°C.

Chromatin Immunoprecipitation Assay (ChIP Assay)

Bejerano Lab

Day 2**Prep:**

- ChIP Buffer III (@ rm. temp)
- incubator/oven @ 65°C
- TSE I (@ rm. temp)
- TSE II (@ rm. temp)
- TE (@ rm. temp)
- Elution Buffer (fresh, @ rm. temp)

Solutions:1.) **TSE I:** 20 mM Tris-HCl (pH 8.1), 150 mM NaCl, 2 mM EDTA, 0.1% SDS, 1% Triton X-100

1 M Tris-HCl, pH 8.1	1 ml
5 M NaCl	1.5 ml
0.5 M EDTA	200 ul
10% SDS	500 ul
10% Triton X-100	5 ml
<u>ddH₂O</u>	<u>41.8 ml</u>
Total	50 ml

2.) **TSE II:** 20 mM Tris-HCl (pH. 8.1), 500 mM NaCl, 2 mM EDTA, 0.1% SDS, 1% Triton X-100

1 M Tris-HCl, pH 8.1	1 ml
5 M NaCl	5 ml
0.5 M EDTA	200 ul
10% SDS	500 ul
10% Triton X-100	5 ml
<u>ddH₂O</u>	<u>38.3 ml</u>
Total	50 ml

3.) **ChIP Buffer III:** 10 mM Tris-HCl (pH 8.1), 0.25 M LiCl, 1 mM EDTA, 1% NP-40, 1% deoxycholate

1 M Tris-HCl, pH 8.1	500 ul
5 M LiCl	2.5 ml
0.5 M EDTA	100 ul
10% NP-40	5 ml
10% deoxycholic acid	5 ml
<u>ddH₂O</u>	<u>36.9 ml</u>
Total	50 ml

4.) **TE:** 10 mM Tris-HCl (pH 8.1), 1 mM EDTA

1 M Tris-HCl, pH 8.1	500 ul
0.5 M EDTA	100 ul
<u>ddH₂O</u>	<u>49.4 ml</u>

Total 50 ml

5.) **Elution Buffer (make fresh):** 1% SDS, 0.1 M NaHCO₃

10% SDS	100 ul	
1 M NaHCO ₃	100 ul	(84 mg sodium bicarbonate in 1 ml ddH ₂ O, fresh)
ddH ₂ O	800 ul	
Total	1 ml	

Directions:

Immunoprecipitate (do all steps @ 4°C)

- 1.) Add 100 ul cold, blocked agarose bead slurry to each chromatin sample (do @4°C).
- 2.) Rotate @ 4°C for 45 min.
- 3.) Spin down beads @3.4 x 1000 rpm (1.1 rcf), 30 sec, rm. temp.
- 4.) Aspirate supernatant.

Wash (do all steps @ rm. temp)

Wash Buffer order: **TSE I, TSE II, CHIP Buffer III, TE**

- 1.) Add 1 ml of appropriate Wash Buffer to agarose beads.
- 2.) Rotate @ rm. temp. for 5 min.
- 3.) Spin down @ 3.4 x 1000 rpm (1.1 rcf), 30 sec, rm. temp.
- 4.) Aspirate supernatant.
- 5.) Repeat steps 1-4 for each Wash Buffer.

Elute/Reverse Crosslinks

- 1.) Add 100 ul **fresh** Elution Buffer.
- 2.) Rotate 15 min. @ rm. temp.
- 3.) Spin down @ 3.4 x 1000 rpm (1.1 rcf), 30 sec, rm. temp.
- 4.) Transfer supernatant to new tube.
- 5.) Repeat steps 1-4 combining eluates into 1 tube.
- 6.) Add 8 ul 5M NaCl.
- per 200 ul eluate (192 mM NaCl)
- 7.) Incubate @65°C overnight.

Chromatin Immunoprecipitation Assay (ChIP Assay)

Bejerano Lab

Day 3

Prep:

- Heat block @ 45°C
- 1 M Tris-HCl, pH 7.5
- 0.5 M EDTA
- 10 mg/ml proteinase K
- Qiagen Qiaquick PCR Purification Kit

Directions:

1.) Remove tube from 65°C oven and add:

- 8 ul 1M Tris-HCl, pH7.5
- 4 ul 0.5 M EDTA
- 1 ul 10 mg/ml proteinase K

2.) Incubate @ 45°C for 1 hour.

3.) Follow directions as written in Qiagen's Qiaquick PCR Purification Kit to purify un-crosslinked chromatin.

- Elute into 50 ul Qiagen Elution Buffer